



## Aspirator bottles, HDPE

### Aspirator bottles, HDPE

Natural, narrow-neck storage bottles with or without scale. Storage vessel with drain made of food-safe polyethylene. With UV protection for long life. Resistant to acids and chemicals. Special colour and printing possible. With metal carrying belt. With threaded plug and wide thread 22mm. With white screw cap Ø 65mm. With 2 O-rings/seals. Other sizes available on request.

Please order stopcocks and adapter for stopcocks with R 1/2" separately.



Code	Description	Packaging
LLG09002602	Aspirator bottles, HDPE, without scale, 10 Litres	1 pz.
LLG09002607	Aspirator bottles, HDPE, with scale, 10 Litres	1 pz.
LLG09002603	Aspirator bottles, HDPE, without scale, 25 Litres	1 pz.
LLG09002608	Aspirator bottles, HDPE, with scale, 25 Litres	1 pz.
LLG09002601	Aspirator bottles, HDPE, without scale, 5 Litres	1 pz.
LLG09002606	Aspirator bottles, HDPE, with scale, 5 Litres	1 pz.
LLG06227536	Adapter for stopcocks with R 1/2", Litres	1 pz.
LLG06253481	Angle valve with sleeve connection 19mm (3/4"), Litres	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

