

**Digital force gauge FK****Kern & Sohn****Digital force gauge FK**

Robust Push/Pull force gauge for simple measurement.

- Turnable display: automatic direction identification
- Secure operability due to ergonomic design
- Peak-hold function to capture peaks (value is "frozen" for approx. 10 seconds) or Track function mode for a continuous measurement indication
- Selectable measuring units: N, lb, kg, oz
- Auto-Power-Off
- Can be mounted on all SAUTER test stands

Specifications

Precision: 0.5 % of Max

Internal measuring frequency: 1000 Hz

Overload protection: 200 % of Max

Thread: M8

Dimensions (W x D x H): 195 x 82 x 35 mm

Net weight: approx. 0.72 kg

Power supply: 230 V, 50/60 Hz



Code	Description	Packaging
LLG06270956	Digital force gauge FK, Type FK 25.	1 pz.
LLG06270023	Digital force gauge FK, Type FK 250.	1 pz.
LLG06270957	Digital force gauge FK, Type FK 50.	1 pz.
LLG06255717	Digital force gauge FK, Type FK 500.	1 pz.
LLG06270959	Digital force gauge FK, Type FK 1K.	1 pz.
LLG06270958	Digital force gauge FK, Type FK 100.	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

