

**Liquid-Sampler, single hand operation, stainless steel****Bürkle****Liquid-Sampler, single hand operation, stainless steel**

Sampler for thin to viscous liquids. Easy use with one-handed operation. With discharge funnel for safe filling of sample bottles. Material V4A (1.4404).



Code	Description	Packaging
LLG06092997	Liquid-Sampler, single hand operation, stainless steel, open with thumb, 100 ml, Length 540 mm, Diam. 32 mm	1 pz.
LLG06072639	Liquid-Sampler, single hand operation, stainless steel, open with thumb, 50 ml, Length 460 mm, Diam. 32 mm	1 pz.
LLG06253634	Liquid-Sampler, single hand operation, stainless steel, place onto bottle, 50 ml, Length 455 mm, Diam. 32 mm	1 pz.
LLG07622377	Liquid-Sampler, single hand operation, stainless steel, place onto bottle, 100 ml, Length 535 mm, Diam. 32 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

