

**Spatula, stainless steel****Bürkle****Spatula, stainless steel**

Sample spatula made of stainless steel V4A (1.4404). The highly stable and solid design of the spatula allow it to pierce solid substances or sacks and BigBags. The spatula is made of premium stainless steel, high-gloss polished and completely welded, therefore completely free of undercuts. It is ideally suited for taking samples of powders, granulates and pasty media, especially in food and pharmaceutical applications. Autoclavable, flame sterilisation.



Code	Description	Packaging
LLG06283287	Spatula, stainless steel, 50 ml, Length 30 cm, Imm. depth 188 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

