



## Laboratory mixer, BagMixer 400Series S

### Laboratory mixer, BagMixer® 400Series S

BagMixer® 400 S / 400 SW are a new generation of lab blenders, silent, easy-to-use and with sturdy design.

With a significantly lower noise level during mixing, it is suitable for all types of laboratory applications. For high bacterial extraction. No cross contamination with sterile homogenizer filter bags.

**BagMixer® 400 S** without viewing window, **BagMixer® 400 SW** with viewing window.

- Silent mixing with less than 48 db
- Homogenization of the samples in 60 sec.
- LiquidSensor detects leaking bags
- 4 adjustable speeds
- Digital operating display
- LightCode system indicates the mixed state by coloured light
- White LED light illuminates the mixing chamber

All units have a 3-year warranty and life guarantee for shock absorbers and window door.

**Starter Pack:** 1 BagMixer® 400 S / SW, 500 BagFilter® filter bags, 500 BagLight® plain bags, 50 BagClip® closing clips, 1 BagOpen® bag opener/holder and 1 BagRack® storage rack.

**Classic Pack:** 1 BagMixer® 400 S / SW, 500 BagFilter® filter bags, 500 BagPage® filter bags, 500 BagLight® plain bags, 50 BagClip® closing clips, 1 BagOpen® bag opener/holder, 1 BagRack® storage rack, 1 PetriPile® Petri dish storage rack, 1 BagPipet® pipette and tips and 1000 BagTips® Regular tips for BagPipet®.

**Premium Pack:** 1 BagMixer® 400 S / SW, 1 DiluFlow® gravimetric dilutor, 1000 BagFilter® filter bags, 1000 BagPage® filter bags, 2000 BagLight® plain bags, 50 BagClip® closing clips, 1 BagOpen® bag opener/holder, 2 BagRack® storage rack, 2 PetriPile® Petri dish storage rack, 1 BagPipet® pipet and tips and 1000 BagTips® Regular



tips for BagPipet®.

Code	Description	Packaging
LLG06261704	Laboratory mixer, BagMixer® 400Series S, Type BagMixer® 400 S, Volume 50 ... 400 ml, Width 450 mm, Depth 290 mm, Height 370 mm	1 pezzo
LLG06261705	Laboratory mixer, BagMixer® 400Series S, Type BagMixer® 400 SW, Volume 50 ... 400 ml, Width 450 mm, Depth 290 mm, Height 370 mm	1 pezzo





CARLO ERBA Reagents S.A.S.  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



*CARLO ERBA  
Reagents operates with  
a Certified Quality  
Management System*

