

**Blotting paper****Hahnemühle****Blotting paper**

Recommended for blotting techniques with gels, like Southern, Northern and Western Blots, even Semidry-Blots of proteins

- Made from ultrapure cotton linters and cellulose, entirely without additives to ensure that no contamination will occur during the transfer steps
- Tested specifically for chromatographic and blotting techniques to ensure high and uniform capillary actions
- High wet strength for safe handling

**Further sizes on request.**



Code	Description	Packaging
LLG04006156	Blotting paper, medium absorbency, Width 600 mm, Length 580 mm, Thickness 0.9 mm, Weight 320 g / m <sup>2</sup>	50 pz.
LLG04006157	Blotting paper, high absorbency, Width 600 mm, Length 580 mm, Thickness 1.5 mm, Weight 570 g / m <sup>2</sup>	25 pz.
LLG04006154	Blotting paper, medium absorbency, Width 570 mm, Length 460 mm, Thickness 0.35 mm, Weight 192 g / m <sup>2</sup>	100 pz.
LLG04006155	Blotting paper, medium absorbency, Width 600 mm, Length 580 mm, Thickness 0.35 mm, Weight 192 g / m <sup>2</sup>	100 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

