

**Blotting paper****Hahnemühle****Blotting paper**

Recommended for blotting techniques with gels, like Southern, Northern and Western Blots, even Semidry-Blots of proteins

- Made from ultrapure cotton linters and cellulose, entirely without additives to ensure that no contamination will occur during the transfer steps
- Tested specifically for chromatographic and blotting techniques to ensure high and uniform capillary actions
- High wet strength for safe handling

Further sizes on request.



Code	Description	Packaging
LLG04006156	Blotting paper, medium absorbency, Width 600 mm, Length 580 mm, Thickness 0.9 mm, Weight 320 g / m ²	50 pz.
LLG04006157	Blotting paper, high absorbency, Width 600 mm, Length 580 mm, Thickness 1.5 mm, Weight 570 g / m ²	25 pz.
LLG04006154	Blotting paper, medium absorbency, Width 570 mm, Length 460 mm, Thickness 0.35 mm, Weight 192 g / m ²	100 pz.
LLG04006155	Blotting paper, medium absorbency, Width 600 mm, Length 580 mm, Thickness 0.35 mm, Weight 192 g / m ²	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

