

**Mini centrifuges LLG-uniCFUGE 2 and LLG-uniCFUGE 2/5****LLG Labware****Mini centrifuges LLG-uniCFUGE 2 and LLG-uniCFUGE 2/5**

Compact and reasonably priced mini centrifuges with two rotor types. Speed can be adjusted to 4000 min⁻¹ or 6000 min⁻¹. Powerful motor for fast starting. Quick stop when the lid is opened.

LLG-uniCFUGE 2 with **Universal rotor** for 6 x 1.5/2.0 ml and 6 x 0.5 ml tubes and 2 times 8 x 0.2 ml PCR-strips.

LLG-uniCFUGE 2/5 with **rotor for 4 x 5 ml tubes**.

Scope of supply: centrifuge, mains cable

Specifications

Max. RCF:

Max. speed:

Motor: 1200 x g and 2000 x g

4000 min⁻¹ and 6000 min⁻¹

brushless DC motor

Max. load: 6 x 0.5/1.5/2 ml and 2 x 8-tube PCR strips (0.2 ml) or 4 x 5 ml tubes

Dimensions (diam. x H): 140 x 120 mm

Weight: 0.7 kg

Power supply: 100 ... 240 V, 50/60 Hz

Warranty: 3 years



Code	Description	Packaging
LLG06263510	Mini centrifuges uniCFUGE 2 and uniCFUGE 2/5, uniCFUGE 2 with universal rotor, Plug type EU	1 pezzo
LLG06263511	Mini centrifuges uniCFUGE 2 and uniCFUGE 2/5, uniCFUGE 2 with universal rotor, Plug type UK	1 pz.
LLG06263515	Mini centrifuges uniCFUGE 2 and uniCFUGE 2/5, uniCFUGE 2/5 with rotor for 5 ml tubes, Plug type EU	1 pz.
LLG06263516	Mini centrifuges uniCFUGE 2 and uniCFUGE 2/5, uniCFUGE 2/5 with rotor for 5 ml tubes, Plug type UK	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

