

**Cut-Protection Gloves uvex phynomic C3****uvex****Cut-Protection Gloves uvex phynomic C3**

- lightweight and sensitive cut protection safety glove
 - very good mechanical abrasion resistance thanks to the damp-resistant aqua-polymer foam coating
 - good grip in dry and slightly damp areas
 - good cut protection (level 3) and high tear resistance
 - highly breathable coating
 - outstanding tactile feel
- Colour: sky blue
EN: 388



Code	Description	Packaging
LLG06285575	Cut-Protection Gloves uvex phynomic C3, Size 6	1 pz
LLG06285576	Cut-Protection Gloves uvex phynomic C3, Size 7	1 pz
LLG06285577	Cut-Protection Gloves uvex phynomic C3, Size 8	1 pz
LLG06285578	Cut-Protection Gloves uvex phynomic C3, Size 9	1 pz
LLG06285579	Cut-Protection Gloves uvex phynomic C3, Size 10	1 pz
LLG06285580	Cut-Protection Gloves uvex phynomic C3, Size 11	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

