

**Rapid incinerator****Gestigkeit****Rapid incinerator**

For individual and series incineration of solids and liquids. Upper platform with 8 openings of 34mm diameter, for porcelain crucibles up to 50mm diameter. Upper platform and housing made of stainless steel. With 120 minute timer with audible signal after timed period and separate continuous operation switch. A safety switch switches the mains supply off / on automatically on opening / closing the incinerator chamber. With 1.5 metre cable and earthed plug.

Type SVR/E

With electronic temperature control which allows stepless heating adjustment from 10 to 100%. Complete for use.

Type SVD 95

With digital temperature adjustment and display. Complete for use.



Code	Description	Packaging
LLG09884140	Rapid incinerator, Type SVR/E, dim. (W x D x H) 450 x 310 x 180 mm, Rating 2500 W, Temp. range 950 max. °C	1 pz.
LLG09884142	Rapid incinerator, Type Heating element EH 400 for SVR/E (spare part), dim. (W x D x H) mm, Rating W, Temp. range max. °C	1 pz.
LLG09884143	Rapid incinerator, Type Heating element EH 95 for SVD 95 (spare part), dim. (W x D x H) mm, Rating W, Temp. range max. °C	1 pz.
LLG09884141	Rapid incinerator, Type SVD 95, dim. (W x D x H) 450 x 310 x 180 mm, Rating 2500 W, Temp. range 950 max. °C	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

