

**Electronic single channel pipette Picus / Picus NxT**

Sartorius Lab Instruments

**Electronic single channel pipette Picus® / Picus® NxT**

Picus® / Picus® NxT are compact and lightweight pipettes. Fully electronic operation assures high accuracy and precision of results. Picus® NxT provides distinct advantages for highly regulated laboratories.

Fully electronic liquid handling in the volume range of 0.2 µl to 10 ml.

- Highest level of ergonomics provided by low weight, light electronic tip ejection and comfortable handle design
- Extensive range of pipetting modes reduces the needed pipetting steps and speeds up work
- Electronic brake and piston control system provide outstanding accuracy and repeatability of pipetting results, independent of the user
- Intuitive user interface in four language options: English, French, German, Russian
- Adjustment wheel offers fast volume setting and menu navigation
- Optiload enables perfect tip sealing for accurate delivery from each channel
- Safe-Cone Filters prevent the risk of contamination cost-effectively
- Microwell plate tracker guides the user to pipette into the correct wells
- Calibration adjustment in 1, 2 or 3 points

**Picus® NxT**

- Certificate of accredited 3-point calibration (per ISO 17025 and ISO 8655)
- User programmable pipetting protocols enable the storage of three frequently needed pipetting workflows
- 2-level password protection for stored programs to prevent unauthorized changes (optional)
- Pipette locking, e.g. in case of contamination, increases lab safety by disabling the pipette from use.
- Service and calibration reminders help the users to remember important service dates.
- Repeated blow-out helps to dispense the last droplets of e.g. viscous liquids



Code	Description	Packaging
LLG04658709	Electronic single channel pipette Picus® / Picus® NxT, Type Picus® NxT, 0.2-10 µl	
LLG04658714	Electronic single channel pipette Picus® / Picus® NxT, Type Picus® NxT, 500-10000 µl	
LLG04658711	Electronic single channel pipette Picus® / Picus® NxT, Type Picus® NxT, 10-300 µl	
LLG04658712	Electronic single channel pipette Picus® / Picus® NxT, Type Picus® NxT, 50-1000 µl	
LLG07657739	Electronic single channel pipette Picus® / Picus® NxT, Type Picus®, 10-300 µl	
LLG07657874	Electronic single channel pipette Picus® / Picus® NxT, Type Picus®, 5-120 µl	
LLG04658713	Electronic single channel pipette Picus® / Picus® NxT, Type Picus® NxT, 100-5000 µl	
LLG06261420	Electronic single channel pipette Picus® / Picus® NxT, Type Picus®, 100-5000 µl	
LLG06266285	Electronic single channel pipette Picus® / Picus® NxT, Type Picus®, 0.2-10 µl	



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

