



## Liquid Dispenser, ESD

### Liquid Dispenser, ESD

Ideal for a wide range of applications in chemistry and medical laboratories. Made of antistatic HDPE with Teflon tube inside, for a great chemical resistance, ESD-safe.

With hinged lid. Just touch the dispensing dish with a brush, cloth or other applicator and solvent is brought to the surface the instant it is needed. Stainless steel valves seal contents in bottle until used, preventing spills, dripping and escaping fumes from endangering the user. Fluid pumped into the dish does not drip back into the bottle, keeping stored liquids uncontaminated.



Code	Description	Packaging
LLG06272060	Liquid Dispenser, ESD, 180 ml, Anti-splash fluid dispenser	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

