

**Sticky Swabs, ESD****Ideal-tek****Sticky Swabs, ESD**

Sticky Swabs - ESD safe cleaning stems: Dustfree sticky swabs are specially used and designed for a high degree dust-free environment (class10-100). Anti-static treatment can safely dissipate electrostatic charge. Sticky swabs are applicable for manufacturing and assembling static sensitive products. Property: ESD safe, good absorption, soft and good tenacity, excellent agility, non-volatile, environment-friendly, no thrill odor.

Material:

Tip: PVAC (Polyvinyl acetate), Ethyl acetate - Isopentyl acetate

Handle: PVC



Code	Description	Packaging
LLG06272061	Sticky Swabs, ESD, 1 mm Ø, Cleanroom Class ISO 1-2	25 pz.
LLG06272062	Sticky Swabs, ESD, 2 mm Ø, Cleanroom Class ISO 1-2	25 pz.
LLG06272063	Sticky Swabs, ESD, 3 mm Ø, Cleanroom Class ISO 1-2	25 pz.
LLG06272064	Sticky Swabs, ESD, Kit, box of 75 swabs mixed (25 x 1mm, 25 x 2 mm, 25 x 3 mm), Cleanroom Class ISO 1-2	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

