

**Plastic Probes****Ideal-tek****Plastic Probes**

The carbon fiber-reinforced PEEK tools (polyetheretherketone) are ideal for applications in chemical and biological laboratories. They have a very high chemical resistance and a temperature resistance of up to 300 ° C.

- Stirring rod for the preparation of liquid chemical products
- Spatula for the dosing of chemicals
- Spatula for sampling in microscopy
- Scraper for the scratch-free removal of adhesive residues



Code	Description	Packaging
LLG06272055	Plastic Probes, Curved fine tip and flat strong tip chisel-shape	1 pz.
LLG06272054	Plastic Probes, Fine tip cone-shaped and flat strong tip chisel-shape	1 pz.
LLG06272056	Plastic Probes, Flat fine sharp tip and flat large fine tip chisel-shape both	1 pz.
LLG06272057	Plastic Probes, Kit of 3 plastic probes (1 x 6.272 054, 1 x 6.272 055, 1 x 6.272 056)	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

