

**qPCR-Kits and Regents, Delta-Seek**

The Delta Seek detection kits utilise the sensitivity and speed of hydrolysis probes and qPCR to get the most accurate data as quickly as possible. Each of the over 400 targets are specifically designed to ensure the broadest possible detection profile and detection of all clinically relevant strains and subtypes. All test kits are validated in house on multiple qPCR platforms to ensure cross platform functionality. Kits come with negative and positive control samples.

All Delta Seek reagents are lyophilised simplifying the transportation and storage of the kits. Kits are stable at ambient temperatures for 18 months.

Custom Designs possible upon request.

**Specifications**

**Format: Hydrolysis Probes**

**Dye: FAM (target) and VIC (internal controls), other dye sets are available**

**Shipment: Ambient**

**Shipment state: Lyophilised**





**CARLO ERBA Reagents S.A.S.**  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

