

**Bottle top dispenser VITLAB genius2****VITLAB****Bottle top dispenser VITLAB® genius2**

With recirculation valve and calibration mechanism.

- Quick, easy and precise volume adjustment due to interior scalloped track and practical screw mechanism
- Innovative recirculation valve avoids reagent waste
- No drips after tube has been capped
- Easy calibration adjustment
- High-quality materials (PTFE, PFA, FEP, borosilicate glass, platinum-iridium) ensure good chemical resistance to most acids, solvents and bases
- A positive displacement plunger with PFA seal minimizes the crystallization of liquids and therefore the possibility of the plunger jamming
- Easy to dismantle for cleaning, no tools necessary
- Completely autoclavable at 121 °C
- Can be rotated 360° on bottle
- Certificate of conformity according to DIN 12600 DE-M marked
- Margin of error according to DIN EN ISO 8655-5
- With recirculation valve and calibration mechanism
- With GL 45 thread
- With screwable discharge valve

**Items supplied:** 3 respectively 5 threaded adapters\* made of PP (see chart), 1 telescopic filling tube, 1 mounting tool, operation manual, quality certificate stating all test values



Code	Description	Packaging
LLG06283976	Bottle top dispenser VITLAB® genius2, Volume 0.2 - 2 ml, Grad. 0.05 ml, Accuracy max. vol. 0.5** ( $\pm$ R%), Precision max. vol. 0.1 (£ CV%)	1 pz
LLG06283977	Bottle top dispenser VITLAB® genius2, Volume 0.5 - 5 ml, Grad. 0.10 ml, Accuracy max. vol. 0.5** ( $\pm$ R%), Precision max. vol. 0.1 (£ CV%)	1 pz
LLG06283978	Bottle top dispenser VITLAB® genius2, Volume 1 - 10 ml, Grad. 0.20 ml, Accuracy max. vol. 0.5** ( $\pm$ R%), Precision max. vol. 0.1 (£ CV%)	1 pz
LLG06283979	Bottle top dispenser VITLAB® genius2, Volume 2.5 - 25 ml, Grad. 0.50 ml, Accuracy max. vol. 0.5** ( $\pm$ R%), Precision max. vol. 0.1 (£ CV%)	1 pz
LLG06283980	Bottle top dispenser VITLAB® genius2, Volume 5 - 50 ml, Grad. 1 ml, Accuracy max. vol. 0.5** ( $\pm$ R%), Precision max. vol. 0.1 (£ CV%)	1 pz
LLG06283981	Bottle top dispenser VITLAB® genius2, Volume 10 - 100 ml, Grad. 1 ml, Accuracy max. vol. 0.5** ( $\pm$ R%), Precision max. vol. 0.1 (£ CV%)	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

