



Square Media Bottles Nalgene Typ 2019, PETG, sterile

Square Media Bottles Nalgeneä Typ 2019, PETG, sterile

With white HDPE screw closure. Reduced permeability to CO₂/O₂ which reduces pH-shift. Bottles and closures are radiation-sterilised and pyrogen-free to eliminate costly washing, depyrogenation and autoclaving steps. Heatshrink band around closure and neck provides tamper-evident seal. Packed in shrink-wrapped trays. 2L size has molded-in handgrips and a 53-mm (53B) white closure. Bottles are sterile to 10⁻⁶ SAL, non-pyrogenic, non-cytotoxic and comply with USP Class VI guidelines. Sterile. Transparent. Leakproof. Graduated.



Code	Description	Packaging
LLG09103123	Square Media Bottles Nalgene Typ 2019, PETG, sterile, Type 2019, 125 ml, dim. 110 x 54 x 54 mm, Cover type diam. 38 mm/design 430	24 pz.
LLG09103124	Square Media Bottles Nalgene Typ 2019, PETG, sterile, Type 2019, 250 ml, dim. 146 x 61 x 61 mm, Cover type diam. 38 mm/design 430	24 pz.
LLG09103126	Square Media Bottles Nalgene Typ 2019, PETG, sterile, Type 2019, 1000 ml, dim. 220 x 94 x 94 mm, Cover type diam. 38 mm/design 430	12 pz.
LLG09103121	Square Media Bottles Nalgene Typ 2019, PETG, sterile, Type 2019, 30 ml, dim. 64 x 38 x 38 mm, Cover type diam. 20 mm	24 pz.
LLG09103122	Square Media Bottles Nalgene Typ 2019, PETG, sterile, Type 2019, 60 ml, dim. 82 x 41 x 41 mm, Cover type diam. 24 mm	24 pz.
LLG09103127	Square Media Bottles Nalgene Typ 2019, PETG, sterile, Type 2019, 2000 ml, dim. 271 x 116 x 116 mm, Cover type diam. 53 mm/design 53B	6 pz.
LLG09103125	Square Media Bottles Nalgene Typ 2019, PETG, sterile, Type 2019, 500 ml, dim. 177 x 74 x 74 mm, Cover type diam. 38 mm/design 430	12 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

