

**Sediment Tester Rapido****Gerber****Sediment Tester Rapido**

Flow system Original Gerber with electric vacuum pump and drain hose for rapid suction of the milk samples of 500 ml each. The filter plate is then compared with dirt pictures on a control panel. Especially suitable for serial dirt test determinations. Hinged aluminum bottle, dirt field diam. 18 mm.



Code	Description	Packaging
LLG06284815	Sediment Tester Rapido, Sediment Tester Rapido, 230 VAC	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

