

**Cleanroom Gloves AdvanTech529, nitrile****MAPA****Cleanroom Gloves AdvanTech529, nitrile**

Reinforced mechanical resistance for short-term operations. Manufacturing process guarantees a very low level of contaminants: rates of particles, powder-free, non-volatile residue, extractables; compatible with ISO 5 class (Class M3.5 for International system). 100% nitrile formula enables anti-static properties.

Protection against accidental chemical splashes. Recommended for people who are sensitive to the proteins in natural latex. Thickness: 0.1 mm. Packaging: 100 gloves/bag

-Cat. III

According to EN ISO 374-1:2016/Type B, EN ISO 374-5:2016, EN 421



Code	Description	Packaging
LLG06270501	Cleanroom Gloves AdvanTech529, nitrile, Size 6, Length 300 mm	1 pz.
LLG06270502	Cleanroom Gloves AdvanTech529, nitrile, Size 7, Length 300 mm	1 pz.
LLG06270503	Cleanroom Gloves AdvanTech529, nitrile, Size 8, Length 300 mm	1 pz.
LLG06270505	Cleanroom Gloves AdvanTech529, nitrile, Size 10, Length 300 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

