



Safety Eyeshields uvex super fit CR 9178

Safety Eyeshields uvex super fit CR 9178

- Professional protection of products and humans
- First autoclavable protective spectacle with anti-fog effect
- Sterile eye protection in every situation
- Autoclavable at least 10 times, 121°C, 20 minutes
- Innovative design for maximum functionality and wearer comfort
- Non-allergenic material which is free of harmful substances
- Ultra-thin wraparound lens
- Lightest safety spectacle in the world
- EN 166
- UV 400
- Weight 21 g

coating uvex supravision clean



Code	Description	Packaging
LLG07982028	Safety Eyeshields uvex super fit CR 9178, Colour white/light blue , Lens UV 400, 2C-1,2 W 1 FTKN DIN CE , Coating uvex supravision CR	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

