

**Digital CMOS Microscope Cameras ODC****Kern & Sohn****Digital CMOS Microscope Cameras ODC**

Specialists in microscopy for measurement, counting, documentation, archiving and image processing.

- C-Mount USB-camera with CMOS sensor
- USB-cable, eyepiece adapter and calibration slide included
- Kern-Microscope-VIS Software in 6 languages inclusive, for WINDOWS XP and higher



Code	Description	Packaging
LLG06285908	Digital CMOS Microscope Cameras ODC, Type ODC 832, Resolution 5.1 MP, Sensor format 1/2.5" CMOS, Interface USB 3.0	1 pz.
LLG06285783	Digital CMOS Microscope Cameras ODC, Type ODC 824, Resolution 3.1 MP, Sensor format 1/2" CMOS, Interface USB 2.0	1 pz.
LLG06285907	Digital CMOS Microscope Cameras ODC, Type ODC 831, Resolution 3.1 MP, Sensor format 1/3" CMOS, Interface USB 3.0	1 pz.
LLG06285906	Digital CMOS Microscope Cameras ODC, Type ODC 825, Resolution 5.1 MP, Sensor format 1/2.5" CMOS, Interface USB 2.0	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

