

**Digital Microscope Eyepiece Cameras ODC****Kern & Sohn****Digital Microscope Eyepiece Cameras ODC**

Specialists in microscopy for measurement, counting, documentation, archiving and image processing.

- USB-Eyepiece Camera mit CMOS Sensor
- USB-cable, Eyepiece Adapter and calibration slide
- Kern-Microscope-VIS Software in 6 language, for Windows XP and higher



Code	Description	Packaging
LLG06285909	Digital Microscope Eyepiece Cameras ODC, Type ODC 872, Resolution 1.3 MP, Sensor format 1/3" CMOS, Interface USB 2.0	1 pz.
LLG06285910	Digital Microscope Eyepiece Cameras ODC, Type ODC 874, Resolution 3 MP, Sensor format 1/2.7" CMOS, Interface USB 2.0	1 pz.
LLG06285911	Digital Microscope Eyepiece Cameras ODC, Type ODC 881, Resolution 5 MP, Sensor format 1/2.5" CMOS, Interface USB 3.0	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

