

**Caps for wide-mouth reagent bottles, PP (melamine*)****Kautex****Caps for wide-mouth reagent bottles, PP (melamine*)**

Guarantees the functionality and the tightness of the bottle to protect the contents.

Type a) Screw cap, black with 1.5mm PE foam lining. Ensures that the bottle is tightly and safely sealed.

Type b) Tamper-evident cap, blue, with conical seal and tamper-evident ring. Immediately indicates unauthorised opening. Perfect for determining whether the cap has been removed.

Type c) Tamper-evident cap, red, with PE foam lining and tamper-evident ring. Immediately indicates unauthorised opening. Perfect for determining whether the cap has been removed.



Code	Description	Packaging
LLG09072895	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type a, For bottles 100 ml, Colour black, Thread external diam. 32 mm	1 pz.
LLG09072896	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type a, For bottles 250, 350 ml, Colour black, Thread external diam. 45 mm	1 pz.
LLG09072897	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type a, For bottles 500, 750, 1000 ml, Colour black, Thread external diam. 54 mm	1 pz.
LLG09072898	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type a, For bottles 1500 - 4000 ml, Colour black, Thread external diam. 80 mm	1 pz.
LLG09072893	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type b, For bottles 1500 - 4000 ml, Colour blue, Thread external diam. 80 mm	1 pz.
LLG09072904	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type b, For bottles 50 ml, Colour blue, Thread external diam. 28 mm	1 pz.
LLG09072899	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type c, For bottles 50 ml, Colour red, Thread external diam. 28 mm	1 pz.
LLG09072900	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type c, For bottles 100 ml, Colour red, Thread external diam. 32 mm	1 pz.
LLG09072901	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type c, For bottles 250, 350 ml, Colour red, Thread external diam. 45 mm	1 pz.
LLG09072902	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type c, For bottles 500, 750, 1000 ml, Colour red, Thread external diam. 54 mm	1 pz.
LLG09072903	Caps for wide-mouth reagent bottles, PP (melamine*), Type Type c, For bottles 1500 - 4000 ml, Colour red, Thread external diam. 80 mm	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

