


**Chromatographic columns, PTFE- or Valve Stopcock, DURAN tubing****Lenz****Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing**

Made of DURAN® tubing. Chromatography columns with sintered frit, porosity 0.

Available in three versions:

- with beaded rim and PTFE stopcock
- with NS socket and PTFE stopcock
- with NS socket and needle-valve stopcock (bore 0 - 2.5 mm)



Code	Description	Packaging
LLG06223574	Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing, socket NS 14/23, 35 ml, Int. diam. 15 mm, Length 200 mm, Stopcock Valve	1 pz.
LLG06202418	Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing, beaded rim, 1000 ml, Int. diam. 40 mm, Length 800 mm, Stopcock PTFE	1 pz.
LLG06202416	Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing, beaded rim, 125 ml, Int. diam. 20 mm, Length 400 mm, Stopcock PTFE	1 pz.
LLG06203961	Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing, socket NS 29/32, 125 ml, Int. diam. 20 mm, Length 400 mm, Stopcock PTFE	1 pz.
LLG09025912	Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing, socket NS 14/23, 15 ml, Int. diam. 10 mm, Length 200 mm, Stopcock PTFE	1 pz.
LLG09025913	Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing, socket NS 14/23, 23 ml, Int. diam. 10 mm, Length 300 mm, Stopcock PTFE	1 pz.
LLG06205017	Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing, beaded rim, 35 ml, Int. diam. 15 mm, Length 200 mm, Stopcock PTFE	1 pz.
LLG09025914	Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing, socket NS 14/23, 35 ml, Int. diam. 15 mm, Length 200 mm, Stopcock PTFE	1 pz.
LLG06202417	Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing, beaded rim, 430 ml, Int. diam. 30 mm, Length 600 mm, Stopcock PTFE	1 pz.
LLG06303297	Chromatographic columns, PTFE- or Valve Stopcock, DURAN® tubing, socket NS 29/32, 430 ml, Int. diam. 30 mm, Length 600 mm, Stopcock PTFE	1 pz.



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

