

**CHROMABOND HR-P****MACHEREY-NAGEL****CHROMABOND® HR-P****Polystyrene-divinylbenzene adsorbent resin for SPE**

Basic material: Highly porous polystyrene-divinylbenzene copolymer specific surface 1200 m²/g.

Particle size 50-100 µm very high binding capacity, up to 30 % of adsorbent weight.

Recommended application: aromatic compounds, phenols from water, nitroaromatics from water, pesticides from water, PAHs from oil



Code	Description	Packaging
LLG09003489	CHROMABOND® HR-P, Type HR-P, 3 ml, 200 mg	30 pz.
LLG09003490	CHROMABOND® HR-P, Type HR-P, 3 ml, 500 mg	30 pz.
LLG04003506	CHROMABOND® HR-P, Type HR-P*, 6 ml, 500 mg	30 pz.
LLG04003509	CHROMABOND® HR-P, Type HR-P, 6 ml, 200 mg	30 pz.
LLG06224844	CHROMABOND® HR-P, Type HR-P, 6 ml, 1000 mg	30 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

