



Refrigerated Circulators with MX Temperature Controller

Refrigerated Circulators with MX Temperature Controller

The circulators are available in 4 bath sizes combined with 4 different temperature controllers. The large, 180° freely rotatable display Swivel 180° allows good visibility and control of the values set from different angles.

- Reservoir sizes ranging from 7 liters to 75 liters
- Pressure and suction pump (selected models)
- LidDockä lid holder
- DuraTopä, precision-molded phenolic bath top resists laboratory chemicals
- WhisperCool® Environmental Control System
- Cool Commandä refrigeration system
- Washable air filter
- Large bath opening and drain
- Easy hook-up of fluid connections
- Intuitive setup and operation, featuring multi-language help and screen prompts
- Simplified setup and maintenance
- Over-temperature protection
- Safety classes, according to DIN 12876-1 Class III standards for use with flammable liquids (selected models)
- Working Temperature: From -30° up to +135°C
- Temperature Stability: $\pm 0.07^{\circ}\text{C}$
- 3.25" LCD Display
- Pressure Pump with one pump speed
- Maximum Pressure: 1.8 psi (0.12 bar)
- External circulation: closed loop
- Safety class I
- Single-point calibration capability



Code	Description	Packaging
LLG04658606	Refrigerated Circulators with MX Temperature Controller, 7 l	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

