



Heated Circulating Baths with Advanced Programmable (AP) Temperature Controller

The circulators are available in 4 bath sizes combined with 4 different temperature controllers. The large, 180° freely rotatable display Swivel 180° allows good visibility and control of the values set from different angles.

- Reservoir sizes ranging from 7 liters to 75 liters
- Pressure and suction pump (selected models)
- LidDockä lid holder
- DuraTopä, precision-molded phenolic bath top resists laboratory chemicals
- WhisperCool® Environmental Control System
- Cool Commandä refrigeration system
- Washable air filter
- Large bath opening and drain
- Easy hook-up of fluid connections
- Intuitive setup and operation, featuring multi-language help and screen prompts
- Simplified setup and maintenance
- Over-temperature protection
- Safety classes, according to DIN 12876-1 Class III standards for use with flammable liquids (selected models)
- Working Temperature: Ambient +10 °C to +200 °C
- Temperature Stability: ± 0.005 °C
- Open mode programming
- Intuitive 4.3" (10.9 cm) SmartTouch display
- 11 languages: French, German, Spanish, Chinese, Arabic, English, Portuguese, Russian, Hindi, Italian, Korean
- Variable speed pump with open- and closed-loop external circulation capability
- Pump pressure max. 3.6 psi (0.25 bar)
- Swivel 180° Rotating Controller
- On-board connectivity: USB-A & B, Ethernet, RS-232 / RS-485, remote on/off and external temperature probe
- Event scheduling (time & date) with real-time clock
- Review temperature trends for up to 10 days



- Multiple selectable "home" screens
- On-screen help
- Automatic and/or user-adjustable performance optimization
- 10-point calibration capability

Code	Description	Packaging
LLG04658626	Heated Circulating Baths with Advanced Programmable (AP) Temperature Controller, Bath capacity 7 L, Suction pump cap. 12.2 L/min bar	1 pz.
LLG04658629	Heated Circulating Baths with Advanced Programmable (AP) Temperature Controller, Bath capacity 15 L, Suction pump cap. 12.2 L/min bar	1 pz.
LLG04658632	Heated Circulating Baths with Advanced Programmable (AP) Temperature Controller, Bath capacity 20 L, Suction pump cap. 12.2 L/min bar	1 pz.



Code	Description	Packaging
LLG04658636	Heated Circulating Baths with Advanced Programmable (AP) Temperature Controller, Bath capacity 28 L, Suction pump cap. 12.2 L/min bar	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

