

**Stirring Drives for Microtiter Plates MIXdrive MTP****43953****Stirring Drives for Microtiter Plates MIXdrive MTP**

Ultra-flat, compact, wear-free, powerful magnetic stirring drive for microtiter plates with 6/12/24/96 positions, external control unit, hermetically sealed, fully encapsulated stainless steel housing, IP 68, water-, dust- and germproof, robust construction, easy to clean, maximum resistance, submersible in water- resp. oil baths, up to + 50 °C. For use in CO₂ incubators/ovens, removable centering pins, ideal for integration into robots. Individual dimensions and high temperature version up to +200 °C on request.

Please order control units separately.

Specification

Motor rating output: 10 W

Permissible ambient temperature: -10 ... 50 °C in air, 50 °C submerged in water

Dimensions (W x D x H): 140 x 102 x 26 mm

Weight: 1.1 kg

IP code: IP 68

Power supply: 48 VDC



Code	Description	Packaging
LLG06286229	Magnetic stirring system MIXdrive 24 MTP external control	1 pz.
LLG06286227	Magnetic stirring system MIXdrive 6 MTP external control	1 pz.
LLG06286228	Magnetic stirring system MIXdrive 12 MTP external control	1 pz.
LLG06286230	Magnetic stirring system MIXdrive 96 MTP external control	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

