

**Control Units for Stirring Drive MIXdrive MTP****43953****Control Units for Stirring Drive MIXdrive MTP****MIXcontrol 20 MTP**

Stainless steel housing, digital speed display, 100-2000 rpm, 4-step power setting

MIXcontrol 20 RS232

with RS232 interface

Specifications

Speed range: 100 ... 2000 min⁻¹

Motor rating output: 10 W

Permissible ambient temperature: 0 ... 40 °C

Permissible humidity: 80 %

Display: digital

Dimensions (W x D x H): 200 x 155 x 38 mm

Weight: 1.3 kg

IP code: IP 20

Power supply: 100 ... 240 V / 50 ... 60 Hz



Code	Description	Packaging
LLG06286231	Control unit MIXcontrol MTP for MIXdrive MTP -Please order drive separately-	1 pz.
LLG06286232	Control unit MIXcontrol MTP RS232 for stirring drive bioMIXdrive - Please order drive separately-	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

