

**Microcarrier Spinner Flask Magna Flex**

DWK Life Sciences Inc.

**Microcarrier Spinner Flask Magna Flexä**

Celstir® and Magna Flexä spinner flasks from WHEATON® are high-quality products designed for long life, high efficiency and low maintenance. Celstir® spinner flasks are the ideal environment for culturing / growing a wide variety of cells, such as insect cells, hybridomas and many other suspension adapted cells. Magna Flexä flasks are suitable for microcarrier cultures. Celstir® and Magna Flexä flasks are produced from autoclavable, Type I borosilicate glass. WHEATON® spinner flasks can be used on all magnetic stirrers.

**Scope of supply:** Flask, stirring unit with lid and lids for side arms.

- Designed for microcarrier cultures with a magnetic stirrer
- Gentle stirring action provided by a bulb-shaped glass impeller
- 1:1 headspace ratio for ideal gas exchange
- Autoclavable borosilicate glass that conforms to USP Type I and ASTM E 438 Type I, Class A requirements



| Code        | Description                                                                                                                   | Packaging |
|-------------|-------------------------------------------------------------------------------------------------------------------------------|-----------|
| LLG04658731 | Microcarrier Spinner Flask Magna Flex, 250 ml, dim. (dia. x H) 85 x 175 mm, Cap size, top 51-400, Cap size, sidearms 33-430   | 1 pz.     |
| LLG04658730 | Microcarrier Spinner Flask Magna Flex, 125 ml, dim. (dia. x H) 65 x 155 mm, Cap size, top 51-400, Cap size, sidearms 33-430   | 1 pz.     |
| LLG04658732 | Microcarrier Spinner Flask Magna Flex, 500 ml, dim. (dia. x H) 110 x 190 mm, Cap size, top 100-400, Cap size, sidearms GL 45  | 1 pz.     |
| LLG04658733 | Microcarrier Spinner Flask Magna Flex, 1000 ml, dim. (dia. x H) 130 x 250 mm, Cap size, top 100-400, Cap size, sidearms GL 45 | 1 pz.     |
| LLG04658734 | Microcarrier Spinner Flask Magna Flex, 6000 ml, dim. (dia. x H) 258 x 404 mm, Cap size, top 100-400, Cap size, sidearms GL 45 | 1 pz.     |
| LLG06239455 | Microcarrier Spinner Flask Magna Flex, 3000 ml, dim. (dia. x H) 178 x 341 mm, Cap size, top 100-400, Cap size, sidearms GL 45 | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

