



Accessories for Spinner Flasks Celstir and Magna Flex

DWK Life Sciences Inc.

Accessories for Spinner Flasks Celstir® and Magna Flexä

Further spare parts on request.



Code	Description	Packaging
LLG04658738	Type Spare impeller with lid, For bottles 1000 ml, For Magna Flex	1 pz.
LLG04658750	Type Stainless steel shaft assembly, For bottles 1000 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658759	Type Spare impeller with lid, For bottles 1000 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658765	Type Spare flask, For bottles 1000 ml, For Celstir® / Magna Flex	1 pz.
LLG04658735	Type Spare impeller with lid, For bottles 125 ml, For Magna Flex	1 pz.
LLG04658744	Type Stainless steel shaft assembly, For bottles 125 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658756	Type Spare impeller with lid, For bottles 125 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658763	Type Spare flask, For bottles 125 ml, For Celstir® / Magna Flex	1 pz.
LLG04658761	Type Spare flask, For bottles 25 ml, For Celstir®	1 pz.
LLG04658774	Type Spare flask, For bottles 25 ml, For Jacketed Celstir®	1 pz.
LLG04658746	Type Stainless steel shaft assembly, For bottles 250 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658757	Type Spare impeller with lid, For bottles 250 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658764	Type Spare flask, For bottles 250 ml, For Celstir® / Magna Flex	1 pz.
LLG04658739	Type Spare impeller with lid, For bottles 3000 ml, For Magna Flex	1 pz.
LLG04658751	Type Stainless steel shaft assembly, For bottles 3000 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658766	Type Spare flask, For bottles 3000 ml, For Celstir® / Magna Flex	1 pz.
LLG06241172	Type Spare impeller with lid, For bottles 3000 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658755	Type Spare impeller with lid, For bottles 50 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658762	Type Spare flask, For bottles 50 ml, For Celstir®	1 pz.
LLG04658775	Type Spare flask, For bottles 50 ml, For Jacketed Celstir®	1 pz.



Code	Description	Packaging
LLG04658737	Type Spare impeller with lid, For bottles 500 ml, For Magna Flex	1 pz.
LLG04658748	Type Stainless steel shaft assembly, For bottles 500 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658758	Type Spare impeller with lid, For bottles 500 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658740	Type Spare impeller with lid, For bottles 6000 ml, For Magna Flex	1 pz.
LLG04658760	Type Spare impeller with lid, For bottles 6000 ml, For Celstir® / Jacketed Celstir®	1 pz.
LLG04658767	Type Spare flask, For bottles 8000 ml, For Celstir®	1 pz.
LLG04658736	Type Spare impeller with lid, For bottles 250 ml, For Magna Flex	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

