



Open Bath Systems, Stainless Steel

PolyScience

Open Bath Systems, Stainless Steel

An economical alternative to integrated heating baths for liquid temperatures up to 150°C. Controller bridge rests securely on the bath, yet is easily removed for tank cleaning. Generous bath opening provides ready access to samples and the included bath cover improves stability.

- Working temperatures from ambient +10° up to 150°C
- Advanced Programmable, Advanced Digital or MX Temperature Controller
- 6, 10, 20 or 28 liter reservoir
- DuraTop® Chemical Resistant Bridge, Swivel 180° Rotating Controller (on Advanced Programmable and Advanced Digital models)
- Lidded opening for optional cooling coil (10, 20 and 28 liter models)
- External circulation capability
- Complies with DIN 12876-1 Class I safety requirements for use with non-flammable liquids
- Incl. bath cover

Working temperature range: +10° to +150 °C (+135 °C, MX type)

Temperature stability: ± 0.01 °C (± 0.07 °C, MX type)

Heating capacity: 2.2 kW (1.1 kW, MX type)

Pump capacity: 16.7 l/min / 0.25 bar (AP, AD type)

11.9 l/min / 0.12 bar (MX type)



Code	Description	Packaging
LLG04658639	Open Bath Systems, Stainless Steel, Type AP, Filling volume 6 L, dim. (W x D x H) 340 x 206 x 378 mm, Bath opening 100 x 110 mm	1 pz.
LLG04658640	Open Bath Systems, Stainless Steel, Type AD, Filling volume 6 L, dim. (W x D x H) 340 x 206 x 378 mm, Bath opening 100 x 110 mm	1 pz.
LLG04658641	Open Bath Systems, Stainless Steel, Type MX, Filling volume 6 L, dim. (W x D x H) 340 x 206 x 406 mm, Bath opening 100 x 110 mm	1 pz.
LLG04658648	Open Bath Systems, Stainless Steel, Type AP, Filling volume 28 L, dim. (W x D x H) 531 x 342 x 429 mm, Bath opening 257 x 214 mm	1 pz.
LLG04658649	Open Bath Systems, Stainless Steel, Type AD, Filling volume 28 L, dim. (W x D x H) 531 x 342 x 429 mm, Bath opening 257 x 214 mm	1 pz.
LLG04658650	Open Bath Systems, Stainless Steel, Type MX, Filling volume 28 L, dim. (W x D x H) 531 x 342 x 457 mm, Bath opening 257 x 214 mm	1 pz.
LLG04658645	Open Bath Systems, Stainless Steel, Type AP, Filling volume 20 L, dim. (W x D x H) 531 x 342 x 378 mm, Bath opening 264 x 228 mm	1 pz.
LLG04658646	Open Bath Systems, Stainless Steel, Type AD, Filling volume 20 L, dim. (W x D x H) 531 x 342 x 378 mm, Bath opening 264 x 228 mm	1 pz.
LLG04658647	Open Bath Systems, Stainless Steel, Type MX, Filling volume 20 L, dim. (W x D x H) 531 x 342 x 406 mm, Bath opening 264 x 228 mm	1 pz.
LLG04658642	Open Bath Systems, Stainless Steel, Type AP, Filling volume 10 L, dim. (W x D x H) 353 x 342 x 378 mm, Bath opening 99 x 255 mm	1 pz.
LLG04658643	Open Bath Systems, Stainless Steel, Type AD, Filling volume 10 L, dim. (W x D x H) 353 x 342 x 378 mm, Bath opening 99 x 255 mm	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

