

**Dewar flasks Nalgene, Type 4150, HDPE****Thermo Scientific****Dewar flasks Nalgeneä, Type 4150, HDPE**

With vented, insulating HDPE lid. PE covered handle. Unbreakable and completely safe for short term storage of ice water, dry ice solvents and liquid nitrogen. Also suitable for use as warming baths. Chemical-resistant, reinforced walls, filled with urethane foam, are temperature resistant from -196 to 100°C. 1, 2 and 4 litre flasks have a carrying handle.



Code	Description	Packaging
LLG09031961	Dewar flasks Nalgene, Type 4150, HDPE, Type 4150, 1 Litres, Top diam. 96 mm, Int. height 195 mm	1 pz.
LLG09031970	Dewar flasks Nalgene, Type 4150, HDPE, Type 4150, 10 Litres, Top diam. 198 mm, Int. height 396 mm	1 pz.
LLG09031962	Dewar flasks Nalgene, Type 4150, HDPE, Type 4150, 2 Litres, Top diam. 122 mm, Int. height 231 mm	1 pz.
LLG09031964	Dewar flasks Nalgene, Type 4150, HDPE, Type 4150, 4 Litres, Top diam. 158 mm, Int. height 259 mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

