

**LLG-Electrode stand for 3 electrodes****LLG Labware****LLG-Electrode stand for 3 electrodes**

With cable duct. For three electrodes, 3 opening with diam. 16 mm. Additional opening for temperature sensor. Three joints for high flexibility and adjustment options.

Base plate: 105 x 125 mm

Base plate weight: 0.7 kg

Arm weight: 0.1 kg

Max. arm height: 320 mm



Code	Description	Packaging
LLG06286248	Electrode stand for 3 electrodes, Type Electrode stand	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

