



## Laboratory markers edding 8014/8015F

### Laboratory markers edding 8014/8015 F

Fine nib for clear labelling on small surfaces. Particularly suited for smooth and transparent surfaces e.g. slides, reaction vessels, labware and labels. Quick-drying, smudge-proof, waterproof, lightfast and low-odour. Tested quality: extreme resistance to abrasion and high temperatures - tested and approved by TÜV Saarland.



Code	Description	Packaging
LLG06286280	Laboratory markers edding 8014/8015 F, Type 8015 F, Colour Black, Bullet Nib 0.75 mm	1 pz.
LLG06286279	Laboratory markers edding 8014/8015 F, Type 8014, Colour Black, Bullet Nib 1 mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

