

**Membrane filter holders, PC**

Sartorius Lab Instruments

**Membrane filter holders, PC**

Reusable filter holder, without filter. Filter has to be selected separately. For particle removal, sterile filtration of aqueous solutions or for venting of sterile vessels. Autoclavable. With Luer connection.



| Code        | Description  | Packaging |
|-------------|--|-----------|
| LLG09049241 | Membrane filter holders, PC, Filter diam. 13 mm, Pore size 0.45 µm, Connection Female Luer Lock   Luer Slip, Colour Clear, Sterile - | 12 pz     |
| LLG09049242 | Membrane filter holders, PC, Filter diam. 25 mm, Pore size 0.45 µm, Connection Female Luer Lock   Luer Slip, Colour Clear, Sterile - | 12 pz     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

