

**Membrane filter holders, PC**

Sartorius Lab Instruments

Membrane filter holders, PC

Reusable filter holder, without filter. Filter has to be selected separately. For particle removal, sterile filtration of aqueous solutions or for venting of sterile vessels. Autoclavable. With Luer connection.



Code	Description	Packaging
LLG09049241	Membrane filter holders, PC, Filter diam. 13 mm, Pore size 0.45 µm, Connection Female Luer Lock Luer Slip, Colour Clear, Sterile -	12 pz.
LLG09049242	Membrane filter holders, PC, Filter diam. 25 mm, Pore size 0.45 µm, Connection Female Luer Lock Luer Slip, Colour Clear, Sterile -	12 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

