



Vacuum or pressure filtration apparatus, Typ 16510, PC

Sartorius Lab Instruments

Vacuum or pressure filtration apparatus, Typ 16510, PC

For vacuum or pressure filtration up to 2 bar.

For ultra or sterile filtration.

With collecting bottle 250 ml.

For filter 47mm diameter filters.



Code	Description	Packaging
LLG09049243	Vacuum or pressure filtration apparatus, Typ 16510, PC, Type PC filtration apparatus, Volume 250 ml, Filter size Ø 47mm mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

