

**Syringe filters Spartan, Regenerated cellulose****Whatman****Syringe filters Spartan®, Regenerated cellulose**

Versatile use - ready-to-use filter units with a hydrophilic, low protein-binding membrane made of regenerated cellulose.

First class chemical resistance against the usual aqueous and organic HPLC solvents

-Spartan® 13 and 30 are tested and certified for UV absorbance substances at wavelengths of 210nm and 254nm with water, methanol and acetonitrile;

-Spartan® 13 has an extremely low hold up volume; 10µl, for optimal filtration

Applications:

Ideal for HPLC sample preparation to achieve reproducible results. Carefully controlled manufacture guarantees lowest traces of UV- absorbing, extractable components.



Code	Description	Packaging
LLG09049943	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 13 mm, Pore size 0.2 µm, Housing LLF/Mini-Tip	100 pz.
LLG09049944	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 13 mm, Pore size 0.45 µm, Housing LLF/Mini-Tip	100 pz.
LLG09049950	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 13 mm, Pore size 0.45 µm, Housing LLF/Mini-Tip	500 pz.
LLG09049951	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 13 mm, Pore size 0.2 µm, Housing LLF/Mini-Tip	500 pz.
LLG09049941	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 30 mm, Pore size 0.2 µm, Housing LLF/LM	100 pz.
LLG09049942	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 30 mm, Pore size 0.45 µm, Housing LLF/LM	100 pz.
LLG09049948	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 13 mm, Pore size 0.2 µm, Housing LLF/LM	100 pz.
LLG09049949	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 13 mm, Pore size 0.45 µm, Housing LLF/LM	100 pz.
LLG09049960	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 30 mm, Pore size 0.45 µm, Housing LLF/LM	50 pz.
LLG09049965	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 30 mm, Pore size 0.2 µm, Housing LLF/LM	500 pz.
LLG09049980	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 13 mm, Pore size 0.2 µm, Housing LLF/LM	500 pz.
LLG09049959	Syringe filters Spartan®, Regenerated cellulose, Type non sterile, Filter diam. 30 mm, Pore size 0.45 µm, Housing LLF/LM	500 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

