



Filter beaker glass, DURAN, heavy wall

Filter beaker glass, DURAN®, heavy wall

- With approximate graduations and spout
- Glass type I / neutral glass as per USP, EP and JP
- Autoclavable



Code	Description	Packaging
LLG09051054	Filter beaker glass, DURAN®, heavy wall, Nominal capacity 1000 ml	1 pz
LLG09051063	Filter beaker glass, DURAN®, heavy wall, Nominal capacity 2000 ml	1 pz
LLG09051068	Filter beaker glass, DURAN®, heavy wall, Nominal capacity 3000 ml	1 pz
LLG09051073	Filter beaker glass, DURAN®, heavy wall, Nominal capacity 5000 ml	1 pz
LLG09051091	Filter beaker glass, DURAN®, heavy wall, Nominal capacity 20000 ml	1 pz
LLG09051088	Filter beaker glass, DURAN®, heavy wall, Nominal capacity 15000 ml	1 pz
LLG09051036	Filter beaker glass, DURAN®, heavy wall, Nominal capacity 250 ml	1 pz
LLG09051086	Filter beaker glass, DURAN®, heavy wall, Nominal capacity 10000 ml	1 pz
LLG09051044	Filter beaker glass, DURAN®, heavy wall, Nominal capacity 500 ml	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

