

**Filter flasks, glass DURAN****DWK Life Sciences****Filter flasks, glass DURAN®**

DIN 12476. ISO 6556. Erlenmeyer shape. Heavy wall, for vacuum use. With ground 17.5/26 socket for vacuum tubing with 15mm to 18mm external diameter. These filtering flasks fulfil the regulations of the "equipment and product safety regulations".



Code	Description	Packaging
LLG09051263	Filter flasks, glass DURAN®, 2000 ml, Diam. 166 mm, Neck diam. 60 mm, Height 255 mm	1 pz
LLG09051254	Filter flasks, glass DURAN®, 1000 ml, Diam. 135 mm, Neck diam. 45 mm, Height 230 mm	1 pz
LLG09051236	Filter flasks, glass DURAN®, 250 ml, Diam. 85 mm, Neck diam. 34 mm, Height 155 mm	1 pz
LLG09051244	Filter flasks, glass DURAN®, 500 ml, Diam. 105 mm, Neck diam. 34 mm, Height 185 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

