

**Filter funnels, glass DURAN****DWK Life Sciences****Filter funnels, glass DURAN®**

Glass Type I / neutral glass as per USP, EP and JP. Autoclavable.



Code	Description	Packaging
LLG09052315	Filter funnels, glass DURAN®, 75 ml, Diam. 45 mm, Stem diam. 10 mm, Height 132 mm, Porosity 5	1 pz
LLG09052341	Filter funnels, glass DURAN®, 1000 ml, Diam. 120 mm, Stem diam. 22 mm, Height 270 mm, Porosity 1	1 pz
LLG09052342	Filter funnels, glass DURAN®, 1000 ml, Diam. 120 mm, Stem diam. 22 mm, Height 270 mm, Porosity 2	1 pz
LLG09052343	Filter funnels, glass DURAN®, 1000 ml, Diam. 120 mm, Stem diam. 22 mm, Height 270 mm, Porosity 3	1 pz
LLG09052344	Filter funnels, glass DURAN®, 1000 ml, Diam. 120 mm, Stem diam. 22 mm, Height 270 mm, Porosity 4	1 pz
LLG09052345	Filter funnels, glass DURAN®, 1000 ml, Diam. 120 mm, Stem diam. 22 mm, Height 270 mm, Porosity 5	1 pz
LLG09052321	Filter funnels, glass DURAN®, 125 ml, Diam. 60 mm, Stem diam. 10 mm, Height 140 mm, Porosity 1	1 pz
LLG09052322	Filter funnels, glass DURAN®, 125 ml, Diam. 60 mm, Stem diam. 10 mm, Height 140 mm, Porosity 2	1 pz
LLG09052323	Filter funnels, glass DURAN®, 125 ml, Diam. 60 mm, Stem diam. 10 mm, Height 140 mm, Porosity 3	1 pz
LLG09052324	Filter funnels, glass DURAN®, 125 ml, Diam. 60 mm, Stem diam. 10 mm, Height 140 mm, Porosity 4	1 pz
LLG09052325	Filter funnels, glass DURAN®, 125 ml, Diam. 60 mm, Stem diam. 10 mm, Height 140 mm, Porosity 5	1 pz
LLG09052361	Filter funnels, glass DURAN®, 4000 ml, Diam. 175 mm, Stem diam. 30 mm, Height 425 mm, Porosity 1	1 pz
LLG09052363	Filter funnels, glass DURAN®, 4000 ml, Diam. 175 mm, Stem diam. 30 mm, Height 425 mm, Porosity 3	1 pz
LLG09052364	Filter funnels, glass DURAN®, 4000 ml, Diam. 175 mm, Stem diam. 30 mm, Height 425 mm, Porosity 4	1 pz
LLG09052362	Filter funnels, glass DURAN®, 4000 ml, Diam. 175 mm, Stem diam. 30 mm, Height 425 mm, Porosity 2	1 pz
LLG09052301	Filter funnels, glass DURAN®, 50 ml, Diam. 35 mm, Stem diam. 10 mm, Height 130 mm, Porosity 1	1 pz
LLG09052302	Filter funnels, glass DURAN®, 50 ml, Diam. 35 mm, Stem diam. 10 mm, Height 130 mm, Porosity 2	1 pz
LLG09052303	Filter funnels, glass DURAN®, 50 ml, Diam. 35 mm, Stem diam. 10 mm, Height 130 mm, Porosity 3	1 pz
LLG09052304	Filter funnels, glass DURAN®, 50 ml, Diam. 35 mm, Stem diam. 10 mm, Height 130 mm, Porosity 4	1 pz
LLG09052305	Filter funnels, glass DURAN®, 50 ml, Diam. 35 mm, Stem diam. 10 mm, Height 130 mm, Porosity 5	1 pz



Code	Description	Packaging
LLG09052331	Filter funnels, glass DURAN®, 500 ml, Diam. 95 mm, Stem diam. 22 mm, Height 240 mm, Porosity 1	1 pz
LLG09052332	Filter funnels, glass DURAN®, 500 ml, Diam. 95 mm, Stem diam. 22 mm, Height 240 mm, Porosity 2	1 pz
LLG09052333	Filter funnels, glass DURAN®, 500 ml, Diam. 95 mm, Stem diam. 22 mm, Height 240 mm, Porosity 3	1 pz
LLG09052334	Filter funnels, glass DURAN®, 500 ml, Diam. 95 mm, Stem diam. 22 mm, Height 240 mm, Porosity 4	1 pz
LLG09052335	Filter funnels, glass DURAN®, 500 ml, Diam. 95 mm, Stem diam. 22 mm, Height 240 mm, Porosity 5	1 pz
LLG09052311	Filter funnels, glass DURAN®, 75 ml, Diam. 45 mm, Stem diam. 10 mm, Height 132 mm, Porosity 1	1 pz
LLG09052312	Filter funnels, glass DURAN®, 75 ml, Diam. 45 mm, Stem diam. 10 mm, Height 132 mm, Porosity 2	1 pz
LLG09052313	Filter funnels, glass DURAN®, 75 ml, Diam. 45 mm, Stem diam. 10 mm, Height 132 mm, Porosity 3	1 pz
LLG09052314	Filter funnels, glass DURAN®, 75 ml, Diam. 45 mm, Stem diam. 10 mm, Height 132 mm, Porosity 4	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

