



ORP electrodes SenTix ORP, analogue

ORP electrodes SenTix® ORP, analogue

With glass shaft, Ag/AgCl reference system, connections for AS/DIN, AS/DIN-3, AS/BNC, S7 plug head.

Specifications

Temperature range: 0 ... 100 °C

Length / Diam.: 120 mm / 12 mm

Electrolyte: KCl 3 mol/l

Sensor: Platinum

Sensor shape: Disc (4 mm)

Diaphragm: Ceramic



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG07050767 | ORP electrodes SenTix® ORP, analogue, Type SenTix® ORP | 1 pz |



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
Reagents operates with
a Certified Quality
Management System

