



Membrane filters, Polyamide

Membrane filters, Polyamide

Type 250. Polyamide, white. For filtration of alkaline solutions and organic solvents.



Code	Description	Packaging
LLG09053222	Membrane filters, Polyamide, Diam. 25 mm, Pore size 0.2 µm	100 pz.
LLG09053227	Membrane filters, Polyamide, Diam. 25 mm, Pore size 0.45 µm	100 pz.
LLG09053223	Membrane filters, Polyamide, Diam. 47 mm, Pore size 0.2 µm	100 pz.
LLG09053242	Membrane filters, Polyamide, Diam. 47 mm, Pore size 0.45 µm	100 pz.
LLG09053224	Membrane filters, Polyamide, Diam. 50 mm, Pore size 0.2 µm	100 pz.
LLG09053243	Membrane filters, Polyamide, Diam. 50 mm, Pore size 0.45 µm	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

