

**LLG-Test tube rack, PP****LLG Labware****LLG-Test tube rack, PP**

Made of PP. Stackable. Supplied in three parts and are easy to assemble. Alphanumeric grid referenced on top tier for easy identification of probes. Can be used in water baths. With large lateral marking area. Autoclavable.



Code	Description	Packaging
LLG06286175	Test tube rack, PP, For tubes diam. 30 mm, Colour white, dim. (W x D x H) 252 x 107 x 71 mm, Array 21	1 pz.
LLG06286176	Test tube rack, PP, For tubes diam. 30 mm, Colour blue, dim. (W x D x H) 252 x 107 x 71 mm, Array 21	1 pz.
LLG06286177	Test tube rack, PP, For tubes diam. 30 mm, Colour red, dim. (W x D x H) 252 x 107 x 71 mm, Array 21	1 PZ
LLG06286172	Test tube rack, PP, For tubes diam. 20 mm, Colour white, dim. (W x D x H) 252 x 108 x 72 mm, Array 40	1 pz.
LLG06286173	Test tube rack, PP, For tubes diam. 20 mm, Colour blue, dim. (W x D x H) 252 x 108 x 72 mm, Array 40	1 pz.
LLG06286174	Test tube rack, PP, For tubes diam. 20 mm, Colour red, dim. (W x D x H) 252 x 108 x 72 mm, Array 40	1 pz
LLG06286169	Test tube rack, PP, For tubes diam. 16 mm, Colour white, dim. (W x D x H) 251 x 106 x 71 mm, Array 60	1 pz.
LLG06286170	Test tube rack, PP, For tubes diam. 16 mm, Colour blue, dim. (W x D x H) 251 x 106 x 71 mm, Array 60	1 pz.
LLG06286171	Test tube rack, PP, For tubes diam. 16 mm, Colour red, dim. (W x D x H) 251 x 106 x 71 mm, Array 60	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

