



Gel blotting paper

Gel blotting paper

GB 003, GB 005. High purity. Very absorbent. Easy to use. Suitable for all protein and nucleic acid blotting techniques.



Code	Description	Packaging
LLG09057079	Gel blotting paper, Type GB 005, Width 200 mm, Length 200 mm, Wall thickness 1.2 mm	25 pz.
LLG09057075	Gel blotting paper, Type GB 005, Width 580 mm, Length 580 mm, Wall thickness 1.2 mm	25 pz.
LLG09057067	Gel blotting paper, Type GB 003, Width 600 mm, Length 300 mm, Wall thickness 0.8 mm	25 pz.
LLG09057077	Gel blotting paper, Type GB 003, Width 200 mm, Length 200 mm, Wall thickness 0.8 mm	100 pz.
LLG09057078	Gel blotting paper, Type GB 003, Width 600 mm, Length 580 mm, Wall thickness 0.8 mm	50 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

