

**Magnetic stirrers Microstir**

DWK Life Sciences Inc.

Magnetic stirrers Microstir®

- Ideal for low speed shear sensitive mixing
- 5 to 200 rpm, single and four place units
- Programmable, remote control and alarm capable

Specifications**Stirring quantity (H₂O): 3 l****Speed range: 5 to 200 rpm****Permissible ambient temperature: 15 to 40 °C****Permissible relative moisture: 80 % to 31 °C / 50 % at 40 °C****Material: painted aluminium****Power supply: 200-240 VAC, 50/60 Hz****IP code. IP 21****Warranty: 2 years**

Code	Description	Packaging
LLG06286826	WHEATON-Micro-Stir®, 230VAC magnetic stirrer with single place, UK cable	1 pz.
LLG06286827	WHEATON-Micro-Stir®, 230VAC magnetic stirrer with 4 places, EU cable	1 pz.
LLG06286828	WHEATON-Micro-Stir®, 230VAC magnetic stirrer with 4 places, UK cable	1 pz.
LLG06286825	WHEATON-Micro-Stir®, 230VAC magnetic stirrer with single place, EU cable	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

