

**Digital Microscope RED-50X****Motic****Digital Microscope RED-50X**

Digital WiFi microscope. Stream your microscope image live on any device that works with Android and iOS and move freely through the lab or the classroom to show and share your image.

- Monocular head with Built-in digital WiFi camera CMOS 1/3" chip
- Sensor resolution: 1.3 MP
- WiFi transmission: 1280 x 1024 pixels
- Widefield eyepiece WF10X/18 mm, with pointer
- Reversed triple revolving nosepiece
- Achromatic objectives 4X, 10X, 40X S, lockable
- Coarse focusing system
- Gliding stage with built-in 0.65 N.A. condenser and "one-touch" sample clips
- Disc diaphragm
- Cordless LED illumination 20 mA, 3.5 V, 70 mW with intensity control
- Power supply 100-240 V (CE)
- Dust cover, charger and batteries included
- Motic Images 3.0 for Windows, OSX and Linux
- Motic 4-dot calibration slide
- Operation system: Windows, OSX, Android and iOS.



Code	Description	Packaging
LLG06285860	Digital Microscope RED-50X, Type RED-50X, Monocular	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

