

**Laboratory Microscope RED 220****Motic****Laboratory Microscope RED 220**

Designed for advanced biology students, teachers, doctors and lab technicians.

- Binocular head Siedentopf type, 30° inclined, 360° rotating
- Widefield eyepieces WF10X/20 mm
- Reversed quadruple revolving nosepiece
- Achromatic Semi-Plan objectives SP 4X, 10X, 40X S, 100X S Oil
- Coaxial coarse and fine focusing system with tension adjustment
- Built-in low position coaxial mechanical stage
- Focusable Abbe condenser N.A. 1.25 with iris diaphragm
- LED illumination 3 W with intensity control
- Power supply 100-240 V (CE)
- Dust cover, immersion oil (5 ml) and cord hanger
- Monocular-/Trinocular equipment available



Code	Description	Packaging
LLG06285859	Laboratory Microscope RED 220, Type RED 220, Binocular	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

