

**Cleanroom marker, edding 8011****Edding****Cleanroom marker, edding 8011**

For permanent marking of special cleanroom and laminated papers in dust-free environments.

- Smudge-proof and waterproof
- Can be removed from most slides using isopropanol
- Ready to use
- Round tip
- Without any added xylene or toluene



Code	Description	Packaging
LLG06286142	Cleanroom marker, edding 8011, Colour blue, Line width 0.6 mm	10 pz.
LLG06284876	Cleanroom marker, edding 8011, Colour black, Line width 0.6 mm	10 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

